Instant, Natural & Quality Results that Last without Invasive Surgery





HansBiomed Europe Ltd. Tuition House 27-37 St George's Road, London, UK, SW19 4EU HansBiomed: www.hansbiomed.com/en MINT: www.mintpdo.com MINT UK Training: www.mintpdo.co.uk

MIN

E-mail: info.europe@hansbiomed.uk Stay connected with MINT PDO: D in O @mintlifteurope

What is MINT[™]?

MINT is a lifting thread (facial tissue fixation thread) made of dissolvable material polydioxanone (PDO) with 360°, 3D helical barbs to pull and fixate skin tissue for wrinkle improvement.

It is a PDO(Polydioxanone) to deliver safe and outstanding results.

Developed under patented moulding pressed technology, the helically positioned 360°, 3D barbs have high tensile and anchoring strength to provide superior results.

Patented Moulding Press Technology 360° Helically Moulded 3D Barbs PDO Thread with 360 degree helical PDO Threads made through press barbs to increase contouring effect moulding technology to produce barbs **Cross Section** on the surface of the threads 3D Moulded barbs provide maximum anchoring and tensile strength Helically moulded 3D barbs fix and support the skin tissue from all directions

Moulding Method



Anchoring Strength

The patented, helical moulded barbs of MINT provides powerful anchoring strength



Tensile Property

MINT has high tensile property due to its uncompromised structural integrity from its unique moulding press technology



Global Patent Registration

MINT is a global brand with patents on manufacturing methods registered in 10 countries. It is a global brand with unique, unrivaled technological capabilities by securing patent registrations on the manufacturing method for each country.

Approval of PDO Thread For Facial Lifting

MINT is a product currently and continuously exported to the United Kingdom, EMEA, North America, South America, Asia, Oceania as a PDO thread for facial lifting.



MINT Achievements

Global lifting certifications and registrations acquired with international countries such as FDA, CE, and CFDA to assure excellent technological capabilities and reliable safety.







China CFDA

Europe CE Certification Acquisition

US FDA Registration

Korea MFDS Approval as facial tissue fixation thread Registration

Timeless Rejuvenation with MINT[™] PDO







EASY EASY Pare MINT

Thread Length	Barbs Length	USP Size	Needle Gauge	Package
424mm	340mm	1-0	19G	2Threads/1Pouch/Box
25mm • • • •	170mm	34mm	170r	nm 25mm ≯∢ →



PRODUCTS

FIX				Thread Length Barbs Length USP Size			Cog Design: Multi-direction	
	FIX	FIX W Type Cannula	150mm	95mm	2-0	18G	20Threads/5Pouches/Box	
2000 m. 			25mm		95mm		30mm ≁	
FIX-MINI	NOT STR						Cog Design: Multi-direction	
	FIX-MINI		Thread Length	Barbs Length	USP Size	Cannula Gauge	Package	
		W Type Canpula	60mm	50mm	2-0	19G	20Threads/5Pouches/Box	
Gi nu Herizan Maren		10mm	50mm					
PETIT							Cog Design: Bi-direction	
	PETIT	W Type Cannula	Thread Length	Barbs Length	USP Size	Cannula Gauge	Package	
inclum. magizities			W Type 1. Cannula	140mm	70mm	2-0	19G	20Threads/5Pouches/Box
		25mm 2	0mm 5mm ⊁∙ ⊧∙	45mm	**	45mm →		
UP		\bigcirc	Thread Length	Barbs Length		Cappula Gauge	Package	
UP	UP	L Type	95mm	86mm	2-0	19G	20Threads/5Pouches/Box	
	Cannula	43mm	9mm	43mm				
			Uni-direction		Multi-direction	T		
ТІР							Cog Design: Bi-direction	
	TR		Thread Length	Barbs Length	USP Size	Cannula Gauge	Package	
		L Type	40mm	34mm	2-0	19G	20Threads/5Pouches/Box	
	enan Tara an Panan	Cannula	17mm 6mm 17r	nm				

PRODUCTS

mono

MINT mono is a smooth thread that stimulates facial tissue to produce collagen, which volumizes the skin.



*Package: 10pcs x 5 sterile bag 50pcs in a pouch

Needle Size (Gauge) Needle Length (mm)		Thread Length (mm)	Thread Thickness (USP)
	40mm	50mm	
25G	50mm	70mm	
	60mm	90mm	
	90mm	150mm	
26G	40mm	50mm	5-0
	50mm	70mm	
	60mm	90mm	
	90mm	150mm	
27G	40mm	50mm	

mono screw

MINT mono screw is a smooth thread that is coiled, which stimulates the surrounding soft tissue to effectively regenerate collagen for skin volumizing.



*Package: 10pcs x 5 sterile bag 50pcs in a pouch

Needle Size (Gauge)	Needle Length (mm)	Thread Length (mm)	Thread Thickness (USP)	
270	50mm	70mm	5.0	
2/0	60mm	90mm	5-0	
	25mm	30mm	6-0	
200	40mm	50mm		
290	50mm	70mm		
	60mm	90mm		
30G	25mm	30mm	7-0	
	40mm	50mm		
	50mm	70mm		

twin

MINT twin consists of two MINT mono threads that are twisted together. Its twist allows greater stimulation than a single smooth thread to synthesize collagen, which in turn results in more volume and skin rejuvenation.

twinscrem

MINT twin screw is comprised of twisting two MINT mono screw, which are already coiled. The two closely coiled threads stimulate more soft tissue, allowing greater collagen production and skin rejuvenation.



Needle Length(mm) ______ HUB

*Package: 10pcs x 5 sterile bag 50pcs in a pouch

_

Needle Size (Gauge)	Needle Length (mm)	Thread Length (mm)	Thread Thickness (USP)	
26G	60mm	90mm	7-0	
	90mm	150mm		
27G	40mm	50mm		
	50mm	70mm	7-0	
	60mm	90mm		



*Package: 10pcs x 5 sterile bag 50pcs in a pouch

Needle Size (Gauge)	Needle Length (mm)	Thread Length (mm)	Thread Thickness (USP)
25G	40mm	50mm	
	50mm	70mm	6.0
	60mm	90mm	0.0
	90mm	150mm	
26G	40mm	50mm	7.0
	50mm	70mm	7-0





WHICH AREAS CAN BE TARGETED?

MINT PDO threads are barbed and highly effective for

- Sagging Forehead and Droopy Eyebrows
- 2 Sagging Eyelid or Cat Eye Look
- 3 Rejuvenating and Lifting Mid-face
- 4 Rhino Nose
- **5** Low or Blunt Nasal Tip
- 6 Softening The Nasolabial Folds
- Sagging Jowls
- Minimising Double Chin or Neck
- Sagging and Wrinkled Neck